

1 eMethods. OCT scanning protocol for HC

Center	Scanner	Details	
Berlin	Spectralis SD-OCT	Peripapillary ring scan	12°; $16 \leq \text{ART} \leq 100$; 1536 A-Scans per B scan; 496 pixels per A scan 3.5 mm; $16 \leq \text{ART} \leq 100$, 768 A scans per B scan; 496 pixels per A scan
		Macular scan	25°x30°; ART=15; 61 vertical B-scans; 768 A-scans per B-scan
Düsseldorf	Spectralis SD-OCT	Peripapillary ring scan	12°; ART=100; 1536 A-Scans per B scan
Barcelona	Spectralis SD-OCT	Peripapillary ring scan	12°; ART=100; 1536 A-Scans per B scan
Prague	Spectralis SD-OCT	Peripapillary ring scan	12°; ART=100; 1536 A-Scans per B scan

2 Abbreviations: ART: Automatic Real Time function; HC: healthy control; OCT: optical coherence
3 tomography.